




REMEEX[®]

FEMALE


The Readjustable Remeex System for SUI

A BETTER CHOICE FOR ISD AND REDOES



- One precise turn makes the difference
- High success rate in sling failures and ISD cases.
- Proven Efficacy with 15 year clinical follow up data.
- Lifetime readjustability

Remeex Female Readjustable System



• **Excellent Dynamic Pressure Transmission System ***

• **Flexibility to increase/decrease urethral support while the patient is standing & actively participating in achieving continence.**

CODE	PRODUCT	DESCRIPTION
SREM-01	REMEEX FEMALE	Complete adjustable system for Female SUI. Includes: sling, threads, varitensor, screw-driver, manipulator, 2 disconnectors, 2 needles PAS21-O and 1 handle MPAS16.
AR-01	REMEEX READJUSTMENT KIT	Manipulator + 2 disconnectors for post-op readjustments
AR-01LP	REMEEX READJUSTMENT KIT LONG	Long Manipulator + 2 disconnectors for post-op readjustments during patient's lifetime
AR-03	THREAD PASSER RETROPUBIC 5,0	1 Needle holder Remeex 5mm diameter, 230cm. length
AR-05	T- HANDLE 5,0	1 Handle in T for Remeex Needle holder

RISK INFORMATION

Rev 2020 09

(*) US Federal Law restricts this device to sale by or on the order of a physician.

Prior to the intervention, visit the FDA Public Health Notification link:

<https://www.fda.gov/MedicalDevices/ProductsandMedicalProcedures/ImplantsandProsthetics/UroGynSurgicalMesh/ucm345221.htm>

If patient is pregnant or has been informed about allergy/sensitivity to polypropylene mesh, do not undergo this surgery. Intervention may also be evaluated by the physician in case the patient is planning future pregnancies. The sling to be implanted is a permanent implant, and some potential complications associated with surgical procedures present risks and may require additional surgery that may or may not correct the complications. Some potential complications include, among others: perforations of vessels/nerves/bladder/urethra/colon/bowel and other surrounding structures, pain, irritation, discomfort, foreign body reaction, hematoma, infection, mesh erosion/extrusion/migration, fistula formation/inflammation, pain during intercourse (dyspareunia), adhesions, fecal incontinence, constipation, urinary incontinence/retention/obstruction, de novo urgency urinary incontinence, de novo prolapse, nerve damage, mesh shortening.

www.neomedic.com



The Continence Company

The information contained in this brochure is to be used only as a guide to assist with product selection. The information contained herein is provided without warranties of any kind, either expressed or implied, and Neomedic International disclaims any and all liabilities for typographical, printing, or production errors or changes affecting the products and/or the specifications contained herein. For complete product details, see Instructions For Use.